

TSID	Name	SPEC	List price (USD)	Category	Binding Capacity
C0069	Mag Beads COOH (300 nm)	5 mL	170	Polymer Magnetic Beads	
C0069	Mag Beads COOH (300 nm)	25 mL	710	Polymer Magnetic Beads	
C0070	Mag Beads Polymer COOH (300 nm)	5 mL	170	Polymer Magnetic Beads	
C0070	Mag Beads Polymer COOH (300 nm)	25 mL	710	Polymer Magnetic Beads	
C0071	Mag Beads COOH (1 μm)	5 mL	170	Polymer Magnetic Beads	
C0071	Mag Beads COOH (1 μm)	25 mL	710	Polymer Magnetic Beads	
C0072	Mag Beads COOH (2 μm)	5 mL	170	Polymer Magnetic Beads	
C0072	Mag Beads COOH (2 μm)	25 mL	710	Polymer Magnetic Beads	
C0073	Mag Beads COOH (2.8 μm)	5 mL	170	Polymer Magnetic Beads	
C0073	Mag Beads COOH (2.8 μm)	25 mL	710	Polymer Magnetic Beads	
C0074	Mag Beads COOH High capacity (2.8 μm)	5 mL	170	Polymer Magnetic Beads	
C0074	Mag Beads COOH High capacity (2.8 μm)	25 mL	710	Polymer Magnetic Beads	
C0075	Mag Beads COOH (3 μm)	5 mL	170	Polymer Magnetic Beads	
C0075	Mag Beads COOH (3 μm)	25 mL	710	Polymer Magnetic Beads	
C0076	Mag Beads COOH (5 μm)	5 mL	170	Polymer Magnetic Beads	
C0076	Mag Beads COOH (5 μm)	25 mL	710	Polymer Magnetic Beads	
C0077	Mag Beads NH2 (2 μm)	5 mL	170	Polymer Magnetic Beads	
C0077	Mag Beads NH2 (2 μm)	25 mL	710	Polymer Magnetic Beads	
C0078	Magrose Beads OH (10-30 μm)	5 mL	110	Agarose Magnetic Beads	
C0078	Magrose Beads OH (10-30 μm)	50 mL	580	Agarose Magnetic Beads	
C0079	Magrose Beads OH (10-30 μm, Ultra-suspension)	5 mL	140	Agarose Magnetic Beads	
C0079	Magrose Beads OH (10-30 μm, Ultra-suspension)	50 mL	1085	Agarose Magnetic Beads	
C0080	Magrose Beads OH	10 mL	155	Agarose Magnetic Beads	
C0080	Magrose Beads OH	100 mL	1100	Agarose Magnetic Beads	
C0081	Magrose Beads OH 6B	10 mL	170	Agarose Magnetic Beads	
C0081	Magrose Beads OH 6B	100 mL	1115	Agarose Magnetic Beads	
C0082	Magrose Beads NH2 (30-150 μm)	5 mL	95	Agarose Magnetic Beads	
C0082	Magrose Beads NH2 (30-150 μm)	50 mL	470	Agarose Magnetic Beads	
C0083	Magrose Beads NH2 (10-30 μm)	5 mL	110	Agarose Magnetic Beads	
C0083	Magrose Beads NH2 (10-30 μm)	50 mL	580	Agarose Magnetic Beads	
C0084	Magrose Beads NH2 (10-30 μm, Ultra-suspension)	5 mL	140	Agarose Magnetic Beads	
C0084	Magrose Beads NH2 (10-30 μm, Ultra-suspension)	50 mL	1085	Agarose Magnetic Beads	
C0085	Magrose Beads COOH (10-30 μm)	5 mL	110	Agarose Magnetic Beads	
C0085	Magrose Beads COOH (10-30 μm)	50 mL	580	Agarose Magnetic Beads	
C0086	Magrose Beads COOH (10-30 μm, Ultra-suspension)	5 mL	140	Agarose Magnetic Beads	
C0086	Magrose Beads COOH (10-30 μm, Ultra-suspension)	50 mL	1085	Agarose Magnetic Beads	
C0087	Magrose Beads COOH (30-150 μm)	5 mL	95	Agarose Magnetic Beads	
C0087	Magrose Beads COOH (30-150 μm)	50 mL	470	Agarose Magnetic Beads	
C0088	Mag Beads Dextran	5 mL	1220	Dextran Magnetic Beads	
C0089	Mag Beads Streptavidin (300 nm), Sterile	1 mL	148	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0089	Mag Beads Streptavidin (300 nm), Sterile	5 mL	568	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0089	Mag Beads Streptavidin (300 nm), Sterile	10 mL	1018	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0090	Mag Beads Streptavidin (1 μm)	1 mL	95	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0090	Mag Beads Streptavidin (1 μm)	5 mL	355	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0090	Mag Beads Streptavidin (1 μm)	10 mL	635	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥450 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0091	Mag Beads Streptavidin (2 μm)	1 mL	95	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥350 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0091	Mag Beads Streptavidin (2 μm)	5 mL	355	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥350 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0091	Mag Beads Streptavidin (2 μm)	10 mL	635	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥350 pmol/mgBotin-兔IgG結合能力：≥15 μg/mg
C0092	Mag Beads Streptavidin (2.8 μm)	1 mL	95	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0092	Mag Beads Streptavidin (2.8 μm)	5 mL	355	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0092	Mag Beads Streptavidin (2.8 μm)	10 mL	635	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0093	Mag Beads Streptavidin (3 μm)	1 mL	95	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥200 pmol/mgBotin-兔IgG結合能力：≥5 μg/mg
C0093	Mag Beads Streptavidin (3 μm)	5 mL	355	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥200 pmol/mgBotin-兔IgG結合能力：≥5 μg/mg
C0093	Mag Beads Streptavidin (3 μm)	10 mL	635	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥200 pmol/mgBotin-兔IgG結合能力：≥5 μg/mg
C0094	Mag Beads Streptavidin (5 μm)	1 mL	95	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0094	Mag Beads Streptavidin (5 μm)	5 mL	355	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0094	Mag Beads Streptavidin (5 μm)	10 mL	635	Streptavidin Magnetic Beads	Botin-核酸探針結合能力（24nt）：≥300 pmol/mgBotin-兔IgG結合能力：≥10 μg/mg
C0095	Magrose Beads Streptavidin (Ultra-suspension)	5 mL	355	Streptavidin Magnetic Beads	Botin-兔IgG結合能力：800 μg/mL
C0095	Magrose Beads Streptavidin (Ultra-suspension)	10 mL	635	Streptavidin Magnetic Beads	Botin-兔IgG結合能力：800 μg/mL
C0097	Mag Beads Epoxy	10 mL	940	Preactivated Magnetic Beads	
C0098	Mag Beads NHS (300 nm)	1 mL	135	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0098	Mag Beads NHS (300 nm)	5 mL	475	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0098	Mag Beads NHS (300 nm)	10 mL	850	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0099	Mag Beads NHS Kit (2 μm)	1 mL	135	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0099	Mag Beads NHS Kit (2 μm)	5 mL	475	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0099	Mag Beads NHS Kit (2 μm)	10 mL	850	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0100	Mag Beads NHS (2 μm)	50 mL	2560	Preactivated Magnetic Beads	≥30 μg 兔 IgG/mg 磁珠
C0101	Magrose Beads NHS Kit (10-30 μm)	1 mL	135	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠

C0101	Magrose Beads NHS Kit (10-30 μm)	5 mL	475	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0101	Magrose Beads NHS Kit (10-30 μm)	10 mL	850	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0102	Magrose Beads NHS (10-30 μm)	50 mL	2560	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0103	Magrose Beads NHS (30-150 μm)	1 mL	135	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0103	Magrose Beads NHS (30-150 μm)	5 mL	475	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0103	Magrose Beads NHS (30-150 μm)	10 mL	850	Preactivated Magnetic Beads	20~30 mg 兔 IgG/mL 磁珠
C0104	Protein A/G Immunoprecipitation Kit	20 T	140	Immunoprecipitation Kits	0.5-0.6 mg hIgG/mL
C0104B	Protein A/G Immunomagnetic Beads	1 mL	58	Immunoprecipitation Magnetic Beads	0.5-0.6 mg hIgG/mL
C0104B	Protein A/G Immunomagnetic Beads	5 mL	248	Immunoprecipitation Magnetic Beads	0.5-0.6 mg hIgG/mL
C0105	Protein A Immunoprecipitation Kit	40 T	140	Immunoprecipitation Kits	≥0.9 mg hIgG/mL
C0105B	Protein A Immunomagnetic Beads	1 mL	58	Immunoprecipitation Magnetic Beads	≥0.9 mg hIgG/mL
C0105B	Protein A Immunomagnetic Beads	5 mL	248	Immunoprecipitation Magnetic Beads	≥0.9 mg hIgG/mL
C0106	Magrose Beads Protein A (30-150 μm)	5 mL	155	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0106	Magrose Beads Protein A (30-150 μm)	25 mL	531	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0106	Magrose Beads Protein A (30-150 μm)	50 mL	960	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0107	Magrose Beads Protein A (10-30 μm)	5 mL	155	Protein Purification Magnetic Beads	40-45 hIgG/mL
C0107	Magrose Beads Protein A (10-30 μm)	25 mL	531	Protein Purification Magnetic Beads	40-45 hIgG/mL
C0107	Magrose Beads Protein A (10-30 μm)	50 mL	960	Protein Purification Magnetic Beads	40-45 hIgG/mL
C0108	Magrose Beads Protein G	5 mL	155	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0108	Magrose Beads Protein G	25 mL	535	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0108	Magrose Beads Protein G	50 mL	960	Protein Purification Magnetic Beads	25-30 hIgG/mL
C0109	Magrose Beads IDA-Nickel	5 mL	45	Protein Purification Magnetic Beads	30-40 mg/mL
C0109	Magrose Beads IDA-Nickel	25 mL	160	Protein Purification Magnetic Beads	30-40 mg/mL
C0109	Magrose Beads IDA-Nickel	50 mL	285	Protein Purification Magnetic Beads	30-40 mg/mL
C0110	Magrose Beads IDA-Cobalt	5 mL	45	Protein Purification Magnetic Beads	20-30 mg/mL
C0110	Magrose Beads IDA-Cobalt	25 mL	160	Protein Purification Magnetic Beads	20-30 mg/mL
C0110	Magrose Beads IDA-Cobalt	50 mL	285	Protein Purification Magnetic Beads	20-30 mg/mL
C0111	Magrose Beads NTA-Nickel (Ultra-suspension)	5 mL	45	Protein Purification Magnetic Beads	40-50 mg/mL
C0111	Magrose Beads NTA-Nickel (Ultra-suspension)	20 mL	160	Protein Purification Magnetic Beads	40-50 mg/mL
C0111	Magrose Beads NTA-Nickel (Ultra-suspension)	40 mL	286	Protein Purification Magnetic Beads	40-50 mg/mL
C0113	Magrose Beads Heparin	5 mL	160	Protein Purification Magnetic Beads	2-3 mg 抗凝血酶 III/mL
C0114	Magrose Beads Strep-Tactin	5 mL	315	Protein Purification Magnetic Beads	~7 mg Strep-tag II 蛋白/mL
C0114	Magrose Beads Strep-Tactin	50 mL	2015	Protein Purification Magnetic Beads	~7 mg Strep-tag II 蛋白/mL
C0115	Magrose Beads DEAE	5 mL	160	Protein Purification Magnetic Beads	≥110 mg BSA/mL
C0116	Protein G Immunoprecipitation Kit	40 T	140	Immunoprecipitation Kits	≥0.85 mg hIgG/mL
C0116B	Protein G Immunomagnetic Beads	1 mL	58	Immunoprecipitation Magnetic Beads	≥0.85 mg hIgG/mL
C0116B	Protein G Immunomagnetic Beads	5 mL	248	Immunoprecipitation Magnetic Beads	≥0.85 mg hIgG/mL
C0117	Protein L Immunoprecipitation Kit	40 T	140	Immunoprecipitation Kits	≥0.5 mg hIgG/mL
C0117B	Protein L Immunomagnetic Beads	1 mL	58	Immunoprecipitation Magnetic Beads	≥0.5 mg hIgG/mL
C0117B	Protein L Immunomagnetic Beads	5 mL	248	Immunoprecipitation Magnetic Beads	≥0.5 mg hIgG/mL
C0118	HA Tag Immunoprecipitation Kit	40 T	270	Immunoprecipitation Kits	≥ 0.4 mg HA 標籤蛋白/mL
C0118B	HA Tag Immunomagnetic Beads	1 mL	180	Immunoprecipitation Magnetic Beads	≥ 0.4 mg HA 標籤蛋白/mL
C0118B	HA Tag Immunomagnetic Beads	5 mL	730	Immunoprecipitation Magnetic Beads	≥ 0.4 mg HA 標籤蛋白/mL
C0118BH	HA Tag Immunomagnetic Beads (300 nm)	1 mL	180	Immunoprecipitation Magnetic Beads	≥ 0.6 mg HA 標籤蛋白/mL 磁珠
C0118BH	HA Tag Immunomagnetic Beads (300 nm)	5 mL	730	Immunoprecipitation Magnetic Beads	≥ 0.6 mg HA 標籤蛋白/mL 磁珠
C0119	MYC Tag Immunoprecipitation Kit	40 T	330	Immunoprecipitation Kits	≥ 0.5 mg MYC 標籤蛋白/mL
C0119B	MYC Tag Immunomagnetic Beads	1 mL	218	Immunoprecipitation Magnetic Beads	≥ 0.5 mg MYC 標籤蛋白/mL
C0119B	MYC Tag Immunomagnetic Beads	5 mL	858	Immunoprecipitation Magnetic Beads	≥ 0.5 mg MYC 標籤蛋白/mL
C0120	Flag Tag Immunoprecipitation Kit	40 T	510	Immunoprecipitation Kits	≥ 0.6 mg Flag 標籤蛋白/mL
C0120B	Flag Tag Immunomagnetic Beads	1 mL	338	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0120B	Flag Tag Immunomagnetic Beads	2 mL	568	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0120B	Flag Tag Immunomagnetic Beads	5 mL	1308	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0120BH	Flag Tag Immunomagnetic Beads (300 nm)	1 mL	338	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0120BH	Flag Tag Immunomagnetic Beads (300 nm)	2 mL	568	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0120BH	Flag Tag Immunomagnetic Beads (300 nm)	5 mL	1308	Immunoprecipitation Magnetic Beads	≥ 0.6 mg Flag 標籤蛋白/mL
C0121	GST Tag Immunoprecipitation Kit	40 T	330	Immunoprecipitation Kits	≥ 0.6 mg GST 標籤蛋白/mL
C0121B	GST Tag Immunomagnetic Beads	1 mL	218	Immunoprecipitation Magnetic Beads	≥ 0.6 mg GST 標籤蛋白/mL
C0121B	GST Tag Immunomagnetic Beads	5 mL	858	Immunoprecipitation Magnetic Beads	≥ 0.6 mg GST 標籤蛋白/mL
C0122	GFP Tag Immunoprecipitation Kit	40 T	415	Immunoprecipitation Kits	≥ 0.8 mg GFP 標籤蛋白/mL
C0122B	GFP Tag Immunomagnetic Beads	1 mL	275	Immunoprecipitation Magnetic Beads	≥ 0.8 mg GFP 標籤蛋白/mL
C0122B	GFP Tag Immunomagnetic Beads	5 mL	1155	Immunoprecipitation Magnetic Beads	≥ 0.8 mg GFP 標籤蛋白/mL
C0122B	GFP Tag Immunomagnetic Beads	10 mL	2076	Immunoprecipitation Magnetic Beads	≥ 0.8 mg GFP 標籤蛋白/mL
C0123	His Tag Immunoprecipitation Kit	40 T	375	Immunoprecipitation Kits	≥ 0.5 mg His 標籤蛋白/mL
C0123B	His Tag Immunomagnetic Beads	1 mL	248	Immunoprecipitation Magnetic Beads	≥ 0.5 mg His 標籤蛋白/mL
C0123B	His Tag Immunomagnetic Beads	5 mL	1518	Immunoprecipitation Magnetic Beads	≥ 0.5 mg His 標籤蛋白/mL
C0140	Flag Tag Magrose Beads	1 mL	121	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥1.0 mg Flag標籤蛋白/mL
C0140	Flag Tag Magrose Beads	5 mL	549	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥1.0 mg Flag標籤蛋白/mL
C0141	HA Tag Magrose Beads	1 mL	121	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥1.0 mg HA標籤蛋白/mL
C0141	HA Tag Magrose Beads	5 mL	549	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥1.0 mg HA標籤蛋白/mL

C0142	GST Tag Magrose Beads	2 mL	161	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	5~10 mg GST 蛋白
C0142	GST Tag Magrose Beads	10 mL	837	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	5~10 mg GST 蛋白
C0142	GST Tag Magrose Beads	25 mL	1348	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	5~10 mg GST 蛋白
C0143	MYC Tag Magrose Beads	1 mL	207	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥0.5 mg MYC標籤蛋白/mL
C0143	MYC Tag Magrose Beads	5 mL	697	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥0.5 mg MYC標籤蛋白/mL
C0144	GFP Tag Magrose Beads	1 mL	Inquiry	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	>1 mg GFP標籤蛋白/mL
C0144	GFP Tag Magrose Beads	5 mL	Inquiry	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	>1 mg GFP標籤蛋白/mL
C0144	GFP Tag Magrose Beads	10 mL	Inquiry	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	>1 mg GFP標籤蛋白/mL
C0145	Oligo (dT) Mag Beads	1 mL	119	Oligo (dT) Magnetic Beads	mRNA 結合量：1-2 µg/mg磁珠
C0145	Oligo (dT) Mag Beads	5 mL	472	Oligo (dT) Magnetic Beads	mRNA 結合量：1-2 µg/mg磁珠
C0145	Oligo (dT) Mag Beads	10 mL	826	Oligo (dT) Magnetic Beads	mRNA 結合量：1-2 µg/mg磁珠
C0146	Magrose Beads Protein A/G	1 mL	87	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥2.0 mg hIgG/mL
C0146	Magrose Beads Protein A/G	5 mL	372	Immunoprecipitation Magnetic Beads/Protein Purification Magnetic Beads	≥2.0 mg hIgG/mL